




(012) 326 7771  
(012) 326 3200 (Fax)

  
**BIOCLIN**™  
SOLUTIONS

## SUCT-003: MANUAL SUCTION APPARATUS

(7B-1)

Simple hand-operated suction device for mucus and phlegm removal. Small, compact with a 200 mL plastic jar. Can be operated with one hand. The device comes with a zipper pouch for travel.



### Technical Specifications:

- **Max Negative Pressure > 40kpa**
- **Max Flow > 20l/Min**
- **Liquid Collection Jar: 200ml**
- **Weight: 202g**

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE

 [sales@bioclin.co.za](mailto:sales@bioclin.co.za)

 [www.bioclin.co.za](http://www.bioclin.co.za)

 Unit 4, 490 Barolong Street,  
Icon Industrial Park, Sunderland Ridge,  
Centurion, 0157

  
**BIOCLIN**™  
SOLUTIONS