




(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

SUCT-006: PORTABLE PHLEGM SUCTION UNIT (7E-G1)

A compact, electrical suction machine with **built-in rechargeable lithium battery**. The 7E-G1 is suitable for ambulances, emergency, home or clinical use. It can be operated in a vehicle with its 12 V DC input. The **built-in battery** provides running time up to 60 minutes on a full charge and makes it very convenient for use where electrical outlets are not within easy reach.

The 7E-G1 can be used for phlegm and mucus removal and as an aspirator in hospital OTs. It features an all-new design, strong but lightweight ABS body and an oil-free, maintenance-free suction pump with 25 LPM capacity and a detachable power cord. Maximum negative pressure can be conveniently controlled from the regulator on the front panel. It comes with a 1000 mL polycarbonate jar with advanced over-flow protection, tubing kit and additional bacterial filters.



Technical Specifications:

- Input Power : AC 220 V @ 50 Hz, 90 VA / DC 12 V
- Battery: 12 V lithium
- Pumping Capacity : 25 LPM
- Maximum Negative Pressure : 80 kPa
- Operating Noise : < 50 dB
- Jar Capacity : 1000 mL
- Operating Time: 60 minutes on battery

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS