




(012) 326 7771  
(012) 326 3200 (Fax)

  
**BIOCLIN**™  
SOLUTIONS

### SUCT-008: PEDAL SUCTION APPARATUS (7B)

This is a non-electrical, foot pedal operated suction machine. Suitable for use in remote locations and in cases of occasional requirement. It features a powder-coated metal frame with carry handle and 1000 mL glass jar.



#### Technical Specifications:

- **Max Negative Pressure > 80kPa**
- **Negative pressure back-reduction while being closed: <0.005MPa/10 Min**
- **Liquid Collection Jar: 100ml**
- **Pedal pressure: ≤245N**
- **Weight: 4.2kg**

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,  
Icon Industrial Park, Sunderland Ridge,  
Centurion, 0157

  
**BIOCLIN**™  
SOLUTIONS