




(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

LITTMANN

INFANT STETHOSCOPE

.The **3M Littmann Classic II Infant Stethoscope** is specially designed for use on the smallest patients. It delivers high acoustic sensitivity, through its dual-sided chestpiece, for both high and low frequency sounds. It is used by medical professionals and students alike to identify, listen to, and study heart, lung, and other body sounds for physical assessment and diagnosis in infant patients.



Feature

- Non-chill bell sleeve for greater patient comfort
- Dual-sided chestpiece with small 2.7cm diaphragm is ideal for infant patient
- Versatile chestpiece design provides high acoustic sensibility with both a floating Diaphragm and open bell
- Non-chill rim and diaphragm provide comfort
- Soft-sealing eartips provide an excellent seal and comfortable fir
- Heartset is easily adjusted for individual fir and comfort. Angled eartubes align with ear canels

Spares not included

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E60E

 sales@bioclin.co.za

 www.bioclin.co.za

 Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS