




(012) 326 7771  
(012) 326 3200 (Fax)



## LAMP-150

Headlight LED 5w incl 2x batt and charger

*(KS) again launch a new and cost-effective medical headlight and binocular loupes which offer a wide choice of options to suit most surgical procedures*



0530-1-001-000



### features

#### Ultra-light & free

*Lightweight(200g), flexible padded headband provides free range of motion*

#### Even-light

*Homogeneous spot provides even light, light spot of luminance uniform and distinct with no light pollution from outside*

#### Adjust & display

*Full digital brightness control and LCD screen display offer obviously power volume condition and brightness step*

#### Dual-battery

*Dual Lithium Polymer battery packs provide constant steady brightness and longevity to the headlight*

### Lamp specifications

Illumination source:	5W LED
Illumination intensity:	10000 lux — 40000 lux
Color temperature:	5000 — 6000K
Spot size	Max D = 70mm Min D = 15mm at D = 320mm
Lamp working life:	50,000 hours continuous burn
Battery	3.7V Lithium Battery
Battery run time	6 hours
Battery charging time	3.6 hours
Input voltage	100 — 240 VAC

### Surgical Applications

- Close examination
- Intensive care
- Recovery rooms
- E.N.T
- General Dentistry
- High Level Dentistry
- Neuro Surgery
- Vascular Surgery
- Transplant Surgery

Also available with loupes [ LAMP-151]

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E60E

✉ sales@bioclin.co.za

🌐 www.bioclin.co.za

📍 Unit 4, 490 Barolong Street,  
Icon Industrial Park, Sunderland Ridge,  
Centurion, 0157

