




(012) 326 7771  
(012) 326 3200 (Fax)

  
**BIOCLIN**™  
**SOLUTIONS**

## BLAD-023

SURGICAL SCALPEL BLADE



### STERILE SURGICAL SCALPEL BLADE SIZE 12 (BOX OF 100)

SCALPEL HANDLE NO 3 TO BE USED WITH BLADES NO 10 TO 19  
SCALPEL HANDLE NO 4 TO BE USED WITH BLADES NO 20 TO 36

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE



[sales@bioclin.co.za](mailto:sales@bioclin.co.za)



[www.bioclin.co.za](http://www.bioclin.co.za)



Unit 4, 490 Barolong Street,  
Icon Industrial Park, Sunderland Ridge,  
Centurion, 0157

  
**BIOCLIN**™  
**SOLUTIONS**