



(012) 326 7771
 (012) 326 3200 (Fax)



EQUI - 273

RADIO FREQUENCY ELECTROSURGICAL
 UNIT FOR MINOR MONOPOLAR SURGERY



SURTRON 50 D

is a radiofrequency electro-surgical unit suitable for minor monopolar surgery.

SURTRON 50 D

Through the functions selection makes pure cut, cut coagulation (CUT 1 - CUT 2) and pure coagulation possible. The one-board microcontroller and digital display guarantees consistent safety and reliability under any conditions.

SURTRON 50 D

Offers outstanding professional surgical results due to its ergonomic and superior safety features.

Standard Accessories

Reusable monopolar handle	/1
Assorted electrodes (10 pz)	/1
Cable for rod neutral electrode	/1
Rod neutral electrode	/1
Foot switch	/1
Power supply cable 2mt	/1

CONTROLS

- Patient/Plate circuit monitoring
- Output power monitoring
- Self check control

SAFETY

EN60601-1
 EN60601-1-2
 EN60601-2-2
 Electrical Class: I CF
 MDD 93/42/EC Class: II b

Dimension

Weight	2,5kg
Size	W190mm x H85mm x D239mm

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E60E

sales@bioclin.co.za

www.bioclin.co.za

Unit 4, 490 Barolong Street,
 Icon Industrial Park, Sunderland Ridge,
 Centurion, 0157

