

Introducing the CTG L8



**BIOCLIN™
SOLUTIONS**

MEDCERT CE 0482
EN ISO 13485

L8 12.1" CTG

Item code: EQUI-322

CTG monitor, SLT L8 colour 12.1" TFT display with 6 parameters(FHR, Toco, FM, Maternal HR, NIBP, Spo2), Built-in battery backup

PRODUCT DETAIL:

L8 12.1" CTG

Lightweight, compact fetal monitor (antepartum, intrapartum +fetal movement) with options of Twins + Material

FEATURES:

- 1. light and compact design, simple to use front panel controls**
- 2.12.1" TFT Colour screen ,folding 90 degree**
- 3. Built in Memory + battery**
- 4. With printer**
- 5.A standard patient event marker and a clinical event marking button to separately mark Clinical events**
- 6.Auto Fetal movement**
- 7. Multi-crystals ,wide beam form, high sensitivity ultrasound transducer ,low ultrasound power**
- 8. AC or LI-battery operated**
- 9. More than 12 hours data storage, then can be played and reprinted**
- 10.Build-in interface to the central nurse station**

PACKAGE DIMENSION:

35 X 32 X 8.5CM

