




(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

FIAD-038

Head immobiliser -Adult



Features

- Large ear hole on both sides to monitor patient ear canal
- Can be attached and secured on scoop stretcher by means of straps
- Velcro adjustable for head and chin straps
- Easy to clean and sterilize for re-use
- Non absorbent plastic coated foam for infection control

Available in child size (fiad-084)

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS