



(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

FIAD - 084

HEAD BLOCK / IMMOBILISER
CHILD



Specifications

High-density foam and plastic moulded Head Blocks for stand alone use, as well as with Spine boards and Scoop Stretchers.

Features:

- Large ear holes in each head supports for monitoring patients ear canals for bleeding / drainage.
- Can be attached & Secured to a scoop stretcher / trauma board by means of straps
- Velcro adjustable Head & Chin Straps
- X-Ray translucent with minimal effect on MRI or CT Scan
- Easy to clean and sterilize for re-use
- Non absorbant plastic coated foam for infection control

Dimensions

Unfolded Size (LxVxH) : 42x26x18cm

Folded Size (LxWxH) : -

Self-weight : 1Kg

Load Bearing : -

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS