



(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

FIAD-125

EVAC chair-1 person carrying



Specifications

- The stretcher is design for assisting patient down stairs
- The design allow it to be use by one person
- 4 wheel to move on level grounds
- Includes 2 belts to secure patient
- Light weight, compact , easy to carry and to clean
- Capacity : 159kg

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E60E



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park,Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS