




(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

FIAD-113

Basket stretcher 2 piece



Feature

- 2 piece design to be compact and easy to transport and store
- 12 mounted handle for ease of carrying
- High density PE construction
- Fire retardant
- include : lifting bridles , carry bag and moulded mattress for patient comfort
- patient restrain strap
- capacity : 272kg

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS