




(012) 326 7771  
(012) 326 3200 (Fax)

  
**BIOCLIN**™  
SOLUTIONS

Furn-006

Commode mobile epoxy



## Specification

Can be placed over the top of the toilet by removing bucket

- Epoxy coated frame
- With footrest
- Removable armrest
- Wheel with breaks
- Capacity : 135kg

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE



[sales@bioclin.co.za](mailto:sales@bioclin.co.za)



[www.bioclin.co.za](http://www.bioclin.co.za)



Unit 4, 490 Barolong Street,  
Icon Industrial Park, Sunderland Ridge,  
Centurion, 0157

  
**BIOCLIN**™  
SOLUTIONS