




(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

FURN-277

Bed side locker –LP110



Features

- Pull out drawer front double as handy flat surface
- Retractable towel rails
- Internal shelf
- Large bottom shelf
- Open bottom shoe shelf

Specifications

- Body material : ABS construction
- Dimensions : 500 x 495 x 870mm
- Top size :450 x 460 x 18mm
- Castors : 55mm

Available in BLUE and GREY

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E60E



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS