



(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

HBME - 013

HEMOTEST DIGITAL HB METER



Specifications

Parameter: Hemoglobin

Principle: Optical reflectance

Sample: Capillary or venous whole blood (15-20 ul)

Strip: H12 Hemoglobin Test Strip

Speed: Result in 10 seconds

Range: 4.0g/dL g/dL

Memory: 250 test results

Calibration: Factory Calibrated according to ICSH

Automatization: Self-checking, auto-judging and displaying malfunction

Power: One CR2032 Battery

Optional : HBME - 014 (Test Strips For Hemotest 50's)

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E60E



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS