





(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

INST-110

Spencer Stitch Scissors
12cm 13.196A



Spencer Stitch Scissors are **primarily used for suture removal**. These scissors have a small hook shaped tip on one blade that slides under sutures to slightly lift them before cutting for removal. This hook also holds the suture so it does not slip off the blade of the scissors prior to cutting.

Features

- Stainless steel
- Hook shaped tip
- 12 cm long

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS