



(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

INST-826

Tuning fork set in zipper pouch



Feature

- High quality polished surgical steel tuning forks set
- Base plate for rinne and weber testing
- 5 base frequencies supplied in set :
 - **128Hz**
 - **256Hz**
 - **512Hz**
 - **1024Hz**
 - **2048Hz**
- 128 , 256 & 512Hz tuning forks can be purchased individually

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E60E



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS