



(012) 326 7771  
(012) 326 3200 (Fax)

**BIOCLIN**™  
**SOLUTIONS**

## LITTMANN

### MASTER CARDIOLOGY STETHOSCOPE

The 3M Littmann Master Cardiology Stethoscope is the ultimate acoustic Stethoscope for accurate auscultation. It is purpose built for Cardiologists, combining superb acoustic performance with outstanding user friendliness. This innovative product has a patented single-sided 'Master' chest piece. This allows monitoring of both high and low frequency sounds without turning the chest piece over.



### Feature.

- Excellent acoustic seal and comfortable fit with 3M snap tight soft-sealing eartips.
- Handcrafted solid stainless steel chestpiece.
- Outstanding acoustics: The highest in Littmann mechanical stethoscope
- Tunable diaphragm: Hear high or low frequency sounds by slightly adjusting pressure on the chestpiece
- Dual-lumen tubing: two sound paths in one eliminates the rubbing noise of traditional double tubes
- Headset is easily adjustable for individual fit and comfort. Angle eartubes with canals

**Available in various colours**

**2160**

**2163**

**2161  
Edition**

**2175  
Brass**

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,  
Icon Industrial Park, Sunderland Ridge,  
Centurion, 0157

**BIOCLIN**™  
**SOLUTIONS**