



(012) 326 7771
 (012) 326 3200 (Fax)

Air Pressure Mattresses

Pressure Relief Bubble Pad System



- Easy to use adjustment settings
- Flame retardant and medical grade bubble pad
- Adjustable hangers for most bed types
- **Maximum weight capacity: 100 Kg**

MATT-018

EXCEL 2000/
 DOMUS1

Pump Dimension: 25x 12.5x 8.5cm
 Weight: 1.4 kg
 Mattress Dimension: 198x86x6.4 cm
 Weight: 2.0 kg




Excel 2000/ DOMUS 1 is ideal for prevention and treatment of patients at low risk of pressure ulcers. This light weight and compact alternating pressure pad is composed of 2.5"-high bubble cells, and it is used as an alternative to heavier products such as gel overlay mattresses. The system is great for short term usage and easy to install. With a pressure adjustable pump and extension flaps on top and bottom sides of the bubble pad to prevent the pad from sliding, DOMUS 1 is designed to be a great transition mattress for homecare.

Adjustable Hanger



Fire retardant PVC




Weight Indicator



For stage 1

Stage I



PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE

✉ sales@bioclin.co.za

🌐 www.bioclin.co.za

📍 Unit 4, 490 Barolong Street,
 Icon Industrial Park, Sunderland Ridge,
 Centurion, 0157