




(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

MEAS-028

Measuring tape pull out type

Note:

The start scale of this tape is 5cm, not 0cm, the reason being that it starts at 5cm is that 5cm covers the length of the plastic part of the calliper. It starts from 5cm in order to provide correct results.

One side in centimeters
one side in inches



Feature

- Tape measuring is easy to use
- **Measures the following:**
 - .Arm
 - .Chest
 - .Waist
 - .Calf
 - .Thigh and hips
- Provide accurate data to help you track your progress accurately.
- Push-button retraction locking feature ensure snug measurement.
- It is to measure by your Self.

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E60E



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS