



(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

OXYG-035

Portable oxygen concentrator 5L



Feature

- PSA technology
- Adjustable flow rate
- De-humidity function
- Timing : LCD show records and the accumulated working time of machine
- Flow rate : 0-5LPM
- Purity : 93%+-3%

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E60E



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS