




(012) 326 7771  
(012) 326 3200 (Fax)

  
**BIOCLIN**™  
SOLUTIONS

## PRES-003

Pressure infusion bag  
1000ml disposable



## Features

- Single patient use effectively reduces the cross contamination
- Gauge safety valve prevents over-inflation (320mmHg pressure relief)
- Colour-coded gauge (0-300mmHg) for easy and accurate pressure monitoring
- The hook design for easy pressure bag hanging during infusion.

## Configuration

- Hook
- Air bladder
- Pressure gauge
- Pressure relief valve
- Inflatable bulb
- Air tube

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,  
Icon Industrial Park, Sunderland Ridge,  
Centurion, 0157

  
**BIOCLIN**™  
SOLUTIONS