



(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

RESU-002

RESUSCITATOR KIT SILICONE
CHILD COMPLETE



**MEETS ISO 9001/13485 AND CE0434 STANDARD.
STERILIZATION: AUTOCLAVABLE;
BY E0 GAS; COLD STERILANT.**

COMPLETE WITH ACCESSORIES IN DUST PROOF CASE

Specifications

- . CHILD TYPE SILICONE RESUSCITATOR WITH VALVE CONNECTOR x 1 PC
- . NO. 3 SILICONE MASK x 1 PC OXYGEN RESERVOIR BAG & CONNECTOR x 1 PC
- . OXYGEN TUBING W/ADAPTER x 1 PC
- . AIR WAY: 60, 70, 80MM 1 PC EACH
- . CARRY BOX x 1 PC

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS