



(012) 326 7771
 (012) 326 3200 (Fax)



STET - 007

DOCTORS BUDGET DUAL HEAD
 STETHOSCOPE



Specifications

- Dual Head chestpiece: for adult, the membrane with diameter
- 47mm is sensitive for good acoustics.
- Aluminium alloy chestpiece, good for acoustics.
- Standard single PVC tubing - color is optional.

AVAILABLE IN THE FOLOWING COLORS :

- BLACK
- GREEN
- RED
- BLUE
- PURPLE
- YELLOW

Device quantity in each box (PCS/BOX)	Carton size(mm)	Quantity In Each	Gross/net Carton(PCS/ctn)	Gross/net weight for each Carton(PCS/ctn)
1	195x405x28	520x405x310	100	14.5/13

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. EGOE

sales@bioclin.co.za
 www.bioclin.co.za

Unit 4, 490 Barolong Street,
 Icon Industrial Park,Sunderland Ridge,
 Centurion, 0157

