




(012) 326 7771  
(012) 326 3200 (Fax)

  
**BIOCLIN**™  
**SOLUTIONS**

## STET-40013

Littmann Spare Parts  
Kit Classic II **INFANT**

Quality and reliability in one convenient package



### Quick overview

- ❖ 3M Littmann stethoscope offers the high quality that medical professionals expect, and each replacement part and accessory meets the same demanding standard. Conveniently packaged in kits that correspond to our stethoscope models, Littmann stethoscope.

### Kit includes:

- One floating diaphragm for the stethoscope chest-piece
- Black rim for assembly of the diaphragm to the chest-piece
- Grey rim for assembly of the diaphragm to chest-piece

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E60E



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,  
Icon Industrial Park, Sunderland Ridge,  
Centurion, 0157

  
**BIOCLIN**™  
**SOLUTIONS**