




(012) 326 7771  
(012) 326 3200 (Fax)

  
**BIOCLIN**™  
**SOLUTIONS**

## SUCT - 006

TRANSPORTER SINGLE BOTTLE  
SUCTION UNIT NEW MODERN  
PLASTIC DESIGN WITH  
BATTERY BACKUP



### Features

- Long life piston cylinder which requires no lubrication or daily maintenance produces no pollution to the environment.
- Provided with three kinds of electrical sources, AC100V-240V 50/60Hz, external DC12V, and internal battery, it can be conveniently used in many places.
- The internal battery is charged with constant voltage and limited current while the unit is working or being attached to AC100V-240V 50/60Hz or external DC 12V power supply. The capacity of the battery is indicated by green, yellow and orange indicators.
- Overflow protection device can safely and efficiently prevent the liquid from entering into pump. Bottles made of plastic material are transparent, durable and convenient for cleaning.
- Body made of plastic material is handy and beautiful. The wall-mounting frame make it widely used in the room, on vehicle and conveniently hung on the side of the wheelchair.

### Specifications

- Max Vacuum: 20.08MPa (600mmHg)
- Flow rate :  $\geq 20$  L/min
- Jar capacity: 1000mL
- Input :110VA
- Dimension : 410mmx 205mm x 420mm 9
- Vacuum range : 0.02MPa - Max vacuum
- Noise : $\leq 65$ d8(A)
- Operating voltage: DC 12V or AC 100V-240V 50/60Hz
- Battery duration on full charge:  $\geq 30$ min
- G.W.IN.W : 6.5/4.5kg

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E80E



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,  
Icon Industrial Park, Sunderland Ridge,  
Centurion, 0157

  
**BIOCLIN**™  
**SOLUTIONS**