




(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

WHEE-001

Rainbow wheelchair 18"



Specification

- One size- 18"
- Full parts and service backup
- One colour- Royal Blue
- Foldable
- Removable armrests
- Removable and height adjustable footrests
- Solid 24" rear wheels
- 8 x 1" solid front castors
- Padded, durable upholstery
- Capacity : 100kg

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS