




(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

W7670-01 / 16 & 10

767 ANERIOD BP WITH CUFF

The versatile 767 Wall, Desk and Mobile Aneroid, with improved ergonomics and a new modern design, is the smart choice for accurate and reliable readings throughout your facility.



Available in:

- W7670-01 (wall mounted)
- W7670-16 (Desktop)
- W7670-10 (mobile)

Features

- The recessed dial face and high-contrast red pointer increase visibility and reduce parallax errors
- The laser-engraved dial face provides unmatched accuracy
- Certified accurate to +/- 3 mm Hg
- Jeweled movement contributes to a longer life
- A new 5-leg stand provides improved stability and manoeuvrability (W7670-10 only)
- Comes with an adult FlexiPort® BP cuff that is antimicrobial treated to help prevent mildew and bacteria growth
- Latex and Mercury Free

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E60E



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS