



(012) 326 7771  
(012) 326 3200 (Fax)

  
**BIOCLIN**™  
SOLUTIONS

# EQUI-084

FINGERTRIP PULSE OXIMETER  
MD300C29



- High brightness OLED displays SpO<sub>2</sub>, PR and Pulse bar
- Low power consumption; real-time battery status indication
- Weak or unstable signal prompt mark helps for more accuracy
- 2pcs AAA-size batteries
- Automatically power off

## Application Range

Adult, adolescent, child

## SpO<sub>2</sub>

Measurement range: 70%-100%

Resolution: 1%

Measurement accuracy: 70%-100%, ±2%; ≤69%: no definition

## Pulse Rate

Measurement range: 30-250bpm

Resolution: 1bpm

Measurement accuracy: 30-99bpm: ±2bpm; 100-250bpm: ±2%

## DisplayType:

OLED Parameters: SpO<sub>2</sub>, PR, Pulse bar

Mode: 6 display modes

## Dimensions and Weight

Dimension: H1.34" x W1.26" x D2.28"

Weight: 0.11lbs (including batteries)

**Power Requirements** Two AAA-size 1.5V batteries

Power consumption: Less than 40mA

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,  
Icon Industrial Park, Sunderland Ridge,  
Centurion, 0157

  
**BIOCLIN**™  
SOLUTIONS