




(012) 326 7771  
(012) 326 3200 (Fax)

  
**BIOCLIN**™  
SOLUTIONS

## FIAD - 119

SPINE BOARD X-RAY TRANSLUCENT  
PLASTIC HALF LENGTH



### Specifications

- Unique design - light weight and easy to use.
- A total of 14 carry handles all around the stretcher allowing it to be versatile.
- The stretcher is made of high density PE material making the structure rugged and wearable
- Supplied with 3 x quick release safety belts.
- Can be used with head immobilizer at the same time.
- Impervious to all bodily fluids
- X-Ray Translucent
- Maximum Capacity : 79Kg

### Dimension

Unfolded Size (LxWxH): 138x35x5cm  
Folded Size (LxWxH) : -  
Self-weight : 4kg  
Load Bearing : 79kg

Packing Dimension  
Quantity : 1  
Size (LxWxH) :140x38x6cm  
G.W : 5kg

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E60E



[sales@bioclin.co.za](mailto:sales@bioclin.co.za)



[www.bioclin.co.za](http://www.bioclin.co.za)



Unit 4, 490 Barolong Street,  
Icon Industrial Park, Sunderland Ridge,  
Centurion, 0157

  
**BIOCLIN**™  
SOLUTIONS