




(012) 326 7771
(012) 326 3200 (Fax)


BIOCLIN™
SOLUTIONS

FURN-121BL

Overbed table ABS Hi-Lo blue



Specifications:

- Maximum height : 1115mm
- Minimum height : 770mm
- Top size : 90mm x 395mm
- Top material : ABS plastic
- Column material : mild steel frame
- Castors : 50mm
- Infection control : (smooth non-absorbent surface)

(Also available in grey)

- **FURN-121G (Overbed table ABS Hi-Lo GREY)**

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E60E



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157


BIOCLIN™
SOLUTIONS