

# OVERBED TABLE



## Available Colours:

- Blue: FURN-121BL
- Grey: FURN-121G

## SPECIFICATIONS

Specifications	
Maximum height	1115mm
Minimum height	770mm
Top size	900mm x 395mm
Top Material	ABS Plastic
Column Material	Mild steel frame
Castors	50mm
Safe working load	25kg
Product Weight	20kg

Features	
Easy Lift Safety, Fully Mechanical Height Adjustment.	•
Locks in Compression no Release Required for Raising	•

**FURN-121BL**

**FURN-121G**

Non-Tilting, Full top



It is a Lightweight overbed table, easy to manoeuvre, reducing the moving and handling risks.

With a slight upwards pressure, the table-top can be raised or lowered which then locks into position, automatically providing a stable surface.

- **Infection Control** (Smooth non-absorbent surfaces)



(012) 326 7771  
(012) 326 3200 (Fax)



sales@bioclin.co.za



Unit 4, 490 Barolong Street,  
Icon Industrial Park, Sunderland Ridge,  
Centurion, 0157



www.bioclin.co.za

