



RT13-2024

Supply Of Respiratory Aids & Anesthetic Accessories

Contract Period : 36 Months

Validity: 1 December 2024- 30 November 2027

WWW.BIOCLIN.CO.ZA



YOUR
LEVEL 1 BEE
PARTNER



SERVING THE
INDUSTRY
SINCE 1999

**PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E & O E
PLEASE ALSO NOTE THAT PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE**



(012) 326 7771
(012) 326 3200 (Fax)



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157



Venturi Oxygen Mask Colour coded

TENDER ITEM NUMBER:

- RT13-02-001: 24% - Paediatric Colour coded venturi: Blue
- RT13-02-002: 28% - Paediatric Colour coded venturi: Yellow
- RT13-02-003: 40% - Paediatric Colour coded venturi: Pink
- RT13-02-004: 60% - Paediatric Colour coded venturi: Red
- RT13-02-005: 28% - Adult Colour coded venturi: Yellow
- RT13-02-006: 40% - Adult Colour coded venturi: Pink
- RT13-02-007: 60% - Adult Colour coded venturi: Red



TENDER ITEM NUMBER: **RT13-04-004**

Adult nebuliser set consisting of an elongated adult mask, medication chamber and oxygen tubing and swivel connector



Oropharyngeal Airway (guedels)

TENDER ITEM NUMBER:


- RT13-05-003: 40-45mm Size 00 (Pink/Cream) 3.0mm
- RT13-05-004: 50-55mm Size 0 (LightBlue/Grey) 3.5mm
- RT13-05-005: 60-65mm Size 1 (Black/Brown) 4.0mm
- RT13-05-006: 70mm Size 2 (White) 4.0mm
- RT13-05-007: 80mm ID:3 (Green) 4.5mm
- RT13-05-008: 90mm Size 4 (Yellow) 4.5mm
- RT-13-05-009: 100mm Size 5 (Red) 5.0mm



Anaesthetic mask with inflatable rim, Anatomically Correct –SINGLE use.

TENDER ITEM NUMBER:

- RT13-08-025: Infant Size 0
- RT13-08-027: Child Size 2



Intubation Stylet

TENDER ITEM NUMBER:

- RT13-40-002: Paediatric OD: 3.3mm Length: 275mm to fit ETT's 4.0 - 5.5mm 10 FR
- RT13-40-003: Adult OD: 4.7mm Total Length: 335 - 345mm to fit ETT's > 6.0mm 12 FR
- RT13-40-004: Adult extra length OD: 4.7mm Total Length: 365mm to fit ETT's > 6.0mm 14 FR



Solid Intubating bougie Single use

TENDER ITEM NUMBER:

- RT13-42-001: Straight Tip Neonatal OD: 1.7mm Total Length: 500mm to fit ETT's 2.5 - 4.0mm 6 FR
- RT13-42-002: Straight Tip Paediatric OD: 3.3mm Total Length: 600 - 700mm to fit ETT's 4.0 - 6.0mm 10FR
- RT13-42-003: Straight Tip Adult OD: 5.0mm Total Length: 600 - 700mm to fit ETT's > 5.0mm 15 FR
- RT13-42-004: Curved Tip Neonatal OD: 1.7mm Total Length: 500mm to fit ETT's 2.5 - 4.0mm 6 FR
- RT13-42-005: Curved Tip Paediatric OD: 3.3mm Total Length: 600 - 700mm to fit ETT's 4.0 - 6.0mm 10FR
- RT13-42-006: Curved Tip Adult OD: 5.0mm Total Length: 600 - 700mm to fit ETT's > 5.0mm 15 FR

**PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E & O E
PLEASE ALSO NOTE THAT PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE**



(012) 326 7771

(012) 326 3200 (Fax)



sales@bioclin.co.za



www.bioclin.co.za



Unit 4, 490 Barolong Street,
Icon Industrial Park, Sunderland Ridge,
Centurion, 0157