




(012) 326 7771  
(012) 326 3200 (Fax)

  
**BIOCLIN**™  
SOLUTIONS

**SUCT-012:  
ELECTRIC SUCTION  
APPARATUS  
(7A-23A)**

A high-capacity, compact suction machine for general surgical use. The 7A-23A features an oil-free, maintenance-free suction pump with 32 LPM capacity. Maximum negative pressure can be conveniently controlled from the regulator on the front panel. The durable, all-ABS aesthetically designed body features a folding handle. It comes with 2 large 2 L polycarbonate jars, a footswitch, a tubing kit and spare bacterial filters.



**Technical Specifications:**

- **Input Power : AC 220 V @ 50 Hz, 180 VA**
- **Pumping Capacity : 20 LPM**
- **Maximum Negative Pressure : 90 kPa**
- **Operating Noise : < 60 dB**
- **Jar Capacity : 2.5 L X 2**

PICTURES ARE FOR ILLUSTRATIVE PURPOSES ONLY AND MAY VARY FROM ACTUAL PRODUCT. E&OE

 [sales@bioclin.co.za](mailto:sales@bioclin.co.za)

 [www.bioclin.co.za](http://www.bioclin.co.za)

 Unit 4, 490 Barolong Street,  
Icon Industrial Park, Sunderland Ridge,  
Centurion, 0157

  
**BIOCLIN**™  
SOLUTIONS